

File Structure and Backup of NIRS/Systemic/MRS Data

L.Beichert

July 30, 2014

During the piglet experiments at the MRI-lab

Sorts of Data Files and Their Origins

- NIRS-files (.tdf): NIRS work station under C://MyDocuments/LWPxxx/
- systemic data files (.m): main Computer under C://temp/LWPxxx/
- 31p MRS files (excel sheets): ask scanning physicist

Structure for Long-Term Saving

- The data of one experiment LWPxxx from the study studyName is to be saved on the main computer into the folder D://data/studyName/LWPxxx/ as follows:
 - the most recent study log (study_log.xlsx) goes into /data/study-name/ (overwrite older version)
 - the NIRS files go into subfolder LWPxxx/nirs/ (use USB-stick to transfer data from NIRS work station to main PC)
 - the systemic files go into subfolder LWPxxx/systemic/
 - the 31p excel sheets go into subfolder LWPxxx/phosphorus/
- Example:
 - data from experiment LWP330 from the argon-dex study goes into folder D://data/argon-dex/LWP330/
 - NIRS files go into data/study-dir/LWP330/nirs

- systemic files go into data/study-dir/LWP330/systemic
- 3lp excel sheet goes into data/study-dir/LWP330/phosphorus

Backup to Nikoleta-Server

- For a security backup the complete data is uploaded to the Nikoleta-Server regularly. The server contains a folder for each study that is supposed to be identical with the study folder on the local computer.
- How to upload data to the Nikoleta server
 - On the main computer start FileZilla FTP Client
 - Go to File → Site Manager, select 'NikoletaServer', enter username and password and click connect
 - Navigate to the study directory both on the local and the remote site
 - In order to upload a new 'LWP'-folder from the local site right-click on it and select 'Upload' in the context menu
 - After uploading the new data disconnect by clicking 'Server' → 'Disconnect'